

For Food Analysis Explore the Advantages of Planar Chromatography

No. 116, March 2016

CAMAG Bibliography Service
Planar Chromatography
Edited by Gertrud Morlock
cbs@camag.com
published by CAMAG Switzerland

IN THIS ISSUE

Procedures, applications

Determination of the hemolytic activity of saponins by an HPTLC-blood gelatin test..... 2–4

The CAMAG TLC-MS Interface 2 in practice 5–7

Identification of herbal slimming drugs and screening for adulteration by HPTLC..... 9–10

Determination of lactose in foodstuff..... 11–12

Quantification of nicotine in liquids for electronic cigarettes 13–15

Products featured in this issue

CAMAG TLC-MS Interface 2..... 7

CAMAG TLC Scanner 4 16

Column: Know CAMAG

The members of the CAMAG laboratory team 8

CAMAG®

CAMAG (Switzerland)
Sonnenmattstrasse 11 • 4132 Muttenz
Phone +41 61 467 3434 • Fax +41 61 461 0702
info@camag.com • www.camag.com

CAMAG Scientific Inc. (USA)
515 Cornelius Harnett Drive
Wilmington, NC 28401
Phone 800 334 3909 • Fax 910 343 1834
tlc@camag.com • www.camagusa.com

Planar Chromatography in Practice

Determination of the hemolytic activity of saponins by an HPTLC-blood gelatin test




Ms. Valeria Reim, M.Sc. Nutritional Science, is a Ph.D. student in the working group of Prof. Dr. Sascha Rohn in the Institute of Food Chemistry at the University of Hamburg. Within the scope of her research project she deals with the characterization of saponins in plant materials with regard to their hemolytic properties using an *in vitro*-HPTLC-blood gelatin test.

Introduction

Saponins are a heterogeneous group of naturally occurring glycosides of triterpenoids, steroids, and steroid alkaloids. They are ubiquitously distributed in the plant kingdom. Legumes, *i.e.* peas and beans, are nutritionally important sources for the uptake of triterpenoid saponins (e.g. soyasaponin I) whereas oats, asparagus, or fenugreek contain steroid saponins [1]. Plants of the nightshade family (*Solanaceae*) such as tomatoes and potatoes contain the so-called solanum alkaloids or steroid alkaloid saponins. Germinated potato sprouts exhibit a particularly high concentration of α -chaconine, one of the most well-known saponins [2]. Due to their amphiphilic characteristics, saponins are able to lyse cell membranes of erythrocytes *in vitro*, therefore being considered as antinutritional factors, so far [3]. For a time, saponin content in plant extracts was semi-quantitatively assayed using the hemolytic features of saponins [1].

In the following application HPTLC was combined with a hemolytic micromethod for the determination of the hemolytic features of saponins following chromatographic separation directly on the HPTLC plate. Besides selected legume saponins (peas and soy bean), glycoalkaloids from potato sprouts were analyzed concerning their potential hemolytic effects towards human blood. Digitonin from the red foxglove *Digitalis purpurea* served as positive control.



Structures of triterpenoid (soyasaponin I), steroid (digitonin) and steroid alkaloid saponins (α -chaconine)

Sample preparation

Steroid alkaloid saponins from potato sprouts (freeze-dried, pulverized) were extracted with methanol in a ratio of 1:10 (w/v) under continuous stirring at room temperature. The supernatant was dried under a stream of nitrogen; the residue was dissolved in bi-distilled water (6 mL) and purified by solid phase extraction [SPE C18, washing with water – methanol (3:2) and elution of saponins with 6 mL methanol]. Subsequently, the extracts were dried under a stream of nitrogen and dissolved in 1 mL methanol. Likewise, legume saponins were extracted from freeze-dried, ground plant material (hulls, saponin mixture from soy beans) with methanol in a ratio of 1:5 (w/v) about 4 h at 50 °C and non-saponin contents were precipitated with ammonium sulfate solution (0.4 M). The extracts were dried under nitrogen flow, the residue was re-suspended in bi-distilled water (6 mL) and further purified by SPE (elution with aqueous methanol solutions with increasing concentrations of 40–100 %, alkalized with 0.017 % aqueous ammonia solution). The saponin-rich fractions (60 % MeOH) were evaporated to dryness and dissolved in 0.25 mL MeOH (60 %, alkalized).

Standard solutions

Ethanolic solutions of digitonin within a concentration range of 0.4–4.0 mg/mL, soyasaponin I ($\geq 95\%$) and a saponin mixture from soy beans at 1 mg/mL, respectively.

Chromatogram layer

HPTLC plates silica gel 60 F₂₅₄ (Merck), 20 x 10 cm

Sample application

Bandwise with Automatic TLC Sampler (ATS 4), max. 20 tracks, band length 7 or 8 mm, distance from left edge 9 or 10 mm, distance from the lower edge 10 mm, application volume between 10 and 20 μ L for sample and standard solutions. Identical application on each half of the plate.

Chromatography

In the Twin Trough Chamber 20 x 10 cm with chamber saturation (without filter paper) for 30 min, with chloroform – methanol – water 6:4:0.9 to a migration distance of 70 mm, air-drying for 1 h at room temperature, cutting the plates into halves with the smartCut

Post-chromatographic derivatization

Visualization of the bands by dipping one half of the plate in *p*-anisaldehyde-sulphuric acid reagent, immersion time 3 s, air-drying for 10 min, heating for 5 min at 70 °C

Blood-gelatin test

Fixation of the zones (2nd plate half) by immersing in 0.5 % polyisobutylmethacrylate solution (in *n*-hexan – chloroform 20:1) for 10 s; coating the plate with 40 mL blood-gelatin [2 % gelatin (w/v) in PBS buffer (pH 7.4) with 3 % human blood preservation] in a petri dish and storage for 12–48 h at 4 °C


Documentation

TLC Visualizer under UV 366 nm and white light

Results and discussion

As saponins do not possess strong chromophores they are only detectable after post-chromatographic derivatization with *p*-anisaldehyde-sulphuric acid reagent under UV 366 nm and white light, cor-

respondingly. Under the influence of ultraviolet light, the steroid alkaloid saponin α -chaconine in potato sprouts occurred purplish-blue and the steroid saponin digitonin appeared yellowish green. With regard to the saponins from pea hulls, the 2,3-dihydro-2,5-dihydroxy-6-methyl-4H-pyran-4-one saponin showed a violet-blue dye whereas soya-saponin I as well as the saponin mixture from soy beans and the soya-saponin I standard were stained in orange. Visualization under white light revealed a brownish colored α -chaconine, while digitonin being racemic occurred as yellow spot. In contrast, saponins from peas and soy beans were all recorded as violet bands.



Chromatogram under UV 366 nm and white light after derivatization with *p*-anisaldehyde-sulphuric acid (D: digitonin, ES1-ES6: pea hulls, KS: potato sprouts, SSI: soya-saponin I standard, S: saponin mixture from soy beans)

One of the crucial biological features of saponins is their ability to lyse cell membranes of erythrocytes, which was detected by a modified microhemolytic method in this case [1]. Exclusively, digitonin and α -chaconine from potato sprouts known to lyse epithelial and red blood cells showed hemolytic zones as expected [4–6]. For a semi-quantitative evaluation of the hemolytic action digitonin was applied in various concentrations. The hemolytic activity of α -chaconine from an extract made of potato sprouts was comparable with the effect of 16–20 μ g digitonin. About 2–5 mg potato alkaloids per kilogram body weight proved to be acutely toxic and the minimum lethal dose amounts are 3–6 mg/kg body weight [2]. Therefore, intoxications to humans resulting from oral administration of glycoalkaloids are rare. In contrast to the studied steroid and steroid alkaloid compounds, the saponins from peas showed no hemolytic effects *in vitro*. When using prolonged incubation of about 48 h, negative results remained, while the bands produced by digitonin became more intense [6]. Hence, relating to the hemolytic aspect of soya-saponins, legumes and legumes-based products can be consumed without hesitation.



Hemolytic zones at white light: sample extracts (top) and digitonin concentration range (bottom, 4–40 μ g/zone; outlier at 20 μ g/zone; 2 zones, because a mixture of mainly digitonin and digaloinin).


Aside from the easy handling, a particular benefit of the presented HPTLC-blood-gelatin screening is the assignment of detected hemolytic zones to single compounds directly on the plate. Known or new saponin compounds can be evaluated concerning their bioactivity. Moreover, a semi-quantitative evaluation of the hemolytic features is ensured using appropriate standard compounds. This technology is extremely useful by giving valuable insights regarding the anti-nutritive status of saponin-rich extracts or foodstuffs, at least *in vitro*.

Further information is available on request from the authors.

- [1] W. A. Oleszek, Saponins In: Natural Food Antimicrobial Systems, CRC Press (2000) 295
- [2] C. Weiss, Ernährungs-Umschau 54 (2007) 474
- [3] G. Francis *et al.*, Brit J Nutr 88 (2002) 587
- [4] E. A. J. Keukens *et al.*, Biochim Biophys Acta 1279 (1996) 243
- [5] E. A. J. Keukens *et al.*, Biochim Biophys Acta 1240 (1995) 216
- [6] O. P. Sharma *et al.*, Food Chem 132 (2012) 671

Contact: Prof. Dr. Sascha Rohn and Valeria Reim, University of Hamburg, Institute of Food Chemistry, Grindelallee 117, 20146 Hamburg, rohn@chemie.uni-hamburg.de, reimv@chemie.uni-hamburg.de

The CAMAG TLC-MS Interface 2 in practice



François Bretin and Dr. Francis Maquin (right)

The Lead Generation to Candidate Realization (LGCR) platform of Sanofi identifies small molecule leads and progresses them up to registration. Dr. Francis Maquin and Mr. François Bretin are members of the R&D department at LGCR Analytical Sciences (AnSci), which is dedicated to supporting their Business Divisions and Therapeutics Units. Having introduced TLC-MS technology in 2009, Sanofi has now acquired one of the first instruments of the second generation in 2015 and has extensively tested the device.

Introduction

TLC-MS is complementary to HPLC-MS and is widely used as a rapid and reliable method for compound confirmation during chemical synthesis.

TLC-MS is especially used when compounds remain on the column due to high polarity or weak solubility or when analyte detection is poor (no chromophore). Thanks to significantly facilitated cleaning features, the TLC-MS Interface 2 can be used for several weeks without any maintenance.

The set-up is augmented by an MS software controlled valve, adding automation to the process and by an Evaporative Light Scattering Detector (ELSD). The latter enables the analyst to distinguish between compounds that were not eluted from the plate from those that were not sufficiently ionized for MS detection.

Chromatogram layer

TLC and HPTLC plates silica gel 60 F₂₅₄ (Merck), 20 × 10 cm. If necessary the plate size was reduced with the smartCut.

Sample preparation

Sample was collected at different stages of the organic synthesis and diluted with an appropriate solvent (e.g. methanol, dimethyl sulfoxide, chloroform).

Sample application

Manually with disposable micropipettes of 5 to 20 μ L volumes

Editor's Note: The Nanomat 4 is recommended to ensure precise positioning without damage to the layer. Especially for polar extracts, small volumes should be applied (<2 μ L to obtain sharp start zones). For a fully-automated and band-wise application the ATS 4 is recommended.

Chromatography


In a twin through chamber, e.g., with mixtures of methanol and dichloromethane – ethyl acetate, or ethyl acetate and heptane – cyclohexane; the ratios depend on the compound mixtures

Derivatization

Compounds with neither UV/Vis-activity nor native fluorescence can be derivatized by non-destructive derivatization reagents, e.g., the primuline or berberine reagents for lipophilic compounds, and directly eluted with the TLC-MS Interface 2 into the MS. However, for destructive derivatizations, e.g., based on strong acidic carbonization reactions, on both plate sides, the outer track was cut, derivatized, and the respective bands were marked by extrapolation.

Recording of MS spectra


The flow rate of the eluent (methanol – water 19:1, 0.35 mL/min) was split by a tee, and 0.15 mL/min were pumped to the MS (Micromass ZQ, Waters) and 0.20 mL/min to the ELSD (Sedere Sedex 85 LT). The TLC-MS Interface 2 was equipped with an oval elution head. The recording of mass spectra was performed in the positive/negative electrospray mode. For evaluation, the MS software (Mass Lynx V4.0/Open Lynx, Waters) was used.



Tee to split the eluent for two detection modes (MS and ELSD)




TLC/HPTLC-ELSD-MS system configuration



Typical chromatogram and TLC-MS spectrum

Results and discussion

From May 2015 to the end of October 2015, the TLC-MS Interface 2 was heavily used without the need for maintenance. We were highly impressed by its performance, as on average 200 elutions were done each month. Compared to the previous version, major improvements were made in cleaning, *i.e.* the elution path is cleaned of matrix particles with compressed air, preventing the system from becoming blocked. An easily accessible filter has been placed in front of the valve, whereas the previous version had a frit inside the elution head. As shown in the image below, the results obtained with the TLC-MS Interface 2 are as satisfactory as with the previous version. We can highly recommend the TLC-MS Interface 2 which allowed us to elute hundreds of zones without any interruption.



A filter is located inline instead of a frit in the elution head

Further information is available on request from the authors.

Contacts:

- Dr. Francis Maquin and Mr. François Bretin, Sanofi-Aventis, Centre de Recherche, 13, quai Jules Guesde, 94403 Vitry-sur-Seine cedex, France, francis.maquin@sanofi-aventis.com
- Mr. Pierre Bernard-Savary, Chromacim SAS, Pommiers la Placette, France
- Mrs. Véronique de Nailly, BCP Instruments, Irigny, France



CAMAG TLC-MS Interface 2

The CAMAG TLC-MS Interface 2 is the second generation of our module for the pioneering concept of hyphenating HPTLC with mass spectrometry. Plate positioning is significantly simplified. The elution head has been modified and an easily accessible, exchangeable filter has been arranged in front of the valve.

Cleaning is facilitated as compared to the previous version, making it highly efficient. By pushing a button, the elution path is cleaned of matrix particles with compressed air, increasing the lifetime of the filter and preventing the system from becoming blocked. Filters are separately available and can be easily replaced without any modification to the elution head.

The CAMAG TLC-MS Interface 2 allows for rapid and contamination-free elution of TLC/HPTLC zones with online transfer to a mass spectrometer. The TLC-MS Interface 2 can be installed plug & play in almost any LC-MS system without adjustments or mass spectrometer modifications. Depending on the MS system, the presence of a substance can be confirmed via its mass spectrum, or for an unknown substance, the respective sum formula can be obtained within a minute.

The members of the CAMAG laboratory team: HPTLC solutions for your analytical problem!



*From left:
Valeria Maire-Widmer (Information Management),
Dr. Eike Reich (Head of Laboratory),
Ilona Trettin (Application Specialist),
Débora Frommenwiler (Scientist, M.Sc.),
Eliezer Ceniviva (Application Specialist);
Tiên Do (Scientific Support Specialist, M.Sc.)*

The members of the CAMAG laboratory team are experienced scientists and technicians dedicated to analytical excellence. They are focused on the advancement of HPTLC as an analytical tool and offer training, application support and feasibility studies to customers worldwide. In many cases HPTLC can be a powerful, pragmatic and cost-efficient solution to an analytical question. We are very proud to introduce below two ladies who have joined the team recently.



Tiên Do joined the team in July of 2014. She holds a master's degree in chemistry (cosmetic formulations, analytics, and quality management) and is currently finishing her PhD thesis ("Evaluation of HPTLC in the analysis of natural extracts") at the University of Nice. She brings experience in applying different chromatographic techniques and problem solving skills to work in interdisciplinary teams. As scientific support specialist Ms. Do provides training and seminars, and develops HPTLC methods in various field of application. Her outgoing personality, generalist's perspective, and ability to think out of the box when approaching an analytical problem are appreciated by customers all over the world.



Ilona Trettin, a chemical-technical and pharmaceutical-technical assistant by training, has practical experience with planar chromatography for more than 20 years, the last 12 of which working in a strictly regulated environment. Joining the team in October of 2014, she complements the existing expertise in the field of HPTLC for quality control of herbal materials. Aside of her responsibility for the experimental section of the new modular in-house HPTLC courses, which are offered several times per year in Muttenz, Ms. Trettin has taken over method development as part of collaborative projects with the European and US Pharmacopoeia.

The abstracts added to this CBS issue cover as usual a broad range of applications. Included also is a thorough review on instrument platforms for TLC/HPTLC (J. Chromatogr. A 1421, 2015, 184–202) presenting a good survey on the current state-of-the-art. Although there is still room for further improvement with regard to miniaturization, the available instrumentation for HPTLC allows a reproducible planar separation either in a standardized or a completely flexible mode of operation.

Some of the included abstracts refer to papers reporting the application of TLC-bioautography. Examples are “A new TLC bioautographic assay for qualitative and quantitative estimation of lipase inhibitors” (Phytochem. Anal. 27, 2015, 5–12) and “The use of the DPPH radical reagent for the detection of radical-scavengers” (J. Planar Chromatogr. 28, 2015, 443–447). Strictly speaking the term bioautography should only be used when biological detection methods are used that involve at least living organelles or cells. The investigation of enzyme inhibitions and the use of the DPPH radical reagent are something of a quite different nature. Even though the current literature does not make a clear distinction between the two approaches and reviewers and editors – as a last instance – do not seem to enforce it. The inclusive term “TLC/HPTLC with effect-directed detection” might be more suitable to describe this broad and still growing field of planar chromatographic research. However be aware that the term “effect-directed analysis” comprises not only effect-directed detections in TLC/HPTLC but also micro-titer plate or cuvette assays.

Dear friends

The current CBS issue puts its focus on food analysis. Anchrom, the Indian distributor of CAMAG products for more than 30 years, made a contribution on herbal slimming drugs in India and their potential adulteration. Other articles deal with the quantifi-




cation of nicotine in liquids for electronic cigarettes, the simple analysis of lactose in food or the determination of the hemolytic activity of saponins in potato sprouts, pea hulls and soy beans by a HPTLC blood gelatin test.

In addition, the R&D Group at Sanofi, after six successful years using the first generation TLC-MS Interface, reports on an application that demonstrates their experience with the new TLC-MS Interface 2.

In every CBS issue we present several applications: Planar Chromatography in Practice. We would very much appreciate feedback from you, dear CBS readers, whether the one or the other of these have prompted you to use this method in your lab or at least consider it. Please communicate your opinion via cbs@camag.com. Thank you in advance.

Insight into all current TLC/HPTLC instrumentation can be accessed by visiting Pittcon 2016, 6–10 March, Atlanta, USA or Analytica 2016, 10–13 May, Munich, Germany or Ilmac 2016, 20–23 September, Basel, Switzerland. Reserve already now the date for the next HPTLC symposium in Berlin on July 5th–7th, 2017.

Kind regards



Gertrud Morlock
cbs@camag.com

THE CBS CLASSIFICATION SYSTEM

- 1. Reviews and books**
 - Books on TLC
 - Books containing one or several chapters on TLC
 - Books containing frequent TLC information spread over several chapters of other information
- 2. Fundamentals, theory and general**
 - General
 - Thermodynamics and theoretical relationship
 - Relationship between structure and chrom. behaviour
 - Measurement of physico-chemical and related values
 - Optimization of solvent systems
 - Validation of methods
- 3. General techniques** (unless they are restricted to the application within one or two classification sections)
 - New apparatus/techniques for sample preparation
 - Separation material
 - New apparatus for sample application/dosage
 - New apparatus/techniques for chromatogram development
 - New apparatus/techniques for pre- or post-chromatographic derivatization
 - New apparatus/techniques for quantitative evaluation
 - New apparatus/techniques for other TLC steps (distinguished from section 4)
- 4. Special techniques**
 - Automation of sample preparation/application
 - Automation of complex chromatogram developing techniques
 - Automation, computer application in quantitative chromatogram evaluation
 - Combination of TLC with other chromatographic techniques
 - Combination of TLC with other (non-chromatographic) techniques...MS, IR...etc.
- 5. Hydrocarbons and halogen derivatives**
 - Aliphatic hydrocarbons
 - Cyclic hydrocarbons
 - Halogen derivatives
 - Complex hydrocarbon mixtures
- 6. Alcohols**
- 7. Phenols**
- 8. Substances containing heterocyclic oxygen**
 - Flavonoids
 - Other compounds with heterocyclic oxygen
- 9. Oxo compounds, ethers and epoxides**
- 10. Carbohydrates**
 - Mono- and oligosaccharides, structural studies
 - Polysaccharides, mucopolysaccharides, lipopolysaccharides
- 11. Organic acids and lipids**
 - Organic acids and simple esters
 - Prostaglandins
 - Lipids and their constituents
 - Lipoproteins and their constituents
 - Glycosphingolipids (gangliosides, sulfatides, neutral glycosphingolipids)
- 12. Organic peroxides**
- 13. Steroids**
 - Pregnane and androstane derivatives
 - Estrogens
 - Sterols
 - Bile acids and alcohols
 - Ecdysones and other insect steroid hormones
- 14. Steroid glycosides, saponins and other terpenoid glycosides**
- 15. Terpenes and other volatile plant ingredients**
 - Terpenes
 - Essential oils
- 16. Nitro and nitroso compounds**
- 17. Amines, amides and related nitrogen compounds**
 - Amines and polyamines
 - Catecholamines and their metabolites
 - Amino derivatives and amides (excluding peptides)
- 18. Amino acids and peptides, chemical structure of proteins**
 - Amino acids and their derivatives
 - Peptides and peptidic proteinous hormones
- 19. Proteins**
- 20. Enzymes**
- 21. Purines, pyrimidines, nucleic acids and their constituents**
 - Purines, pyrimidines, nucleosides, nucleotides
 - Nucleic acids, RNA, DNA
- 22. Alkaloids**
- 23. Other substances containing heterocyclic nitrogen**
 - Porphyryns and other pyrroles
 - Bile pigments
 - Indole derivatives
 - Pyridine derivatives
 - other N-heterocyclic compounds
- 24. Organic sulfur compounds**
- 25. Organic phosphorus compounds** (other than phospholipids)
- 26. Organometallic and related compounds**
 - Organometallic compounds
 - Boranes, silanes and related non-metallic compounds
 - Coordination compounds
- 27. Vitamins and various growth regulators** (non-peptidic)
- 28. Antibiotics, Mycotoxins**
 - Antibiotics
 - Aflatoxins and other mycotoxins
- 29. Pesticides and other agrochemicals**
 - Chlorinated insecticides
 - Phosphorus insecticides
 - Carbamates
 - Herbicides
 - Fungicides
 - Other types of pesticides and various agrochemicals
- 30. Synthetic and natural dyes**
 - Synthetic dyes
 - Chloroplasts and other natural pigments
- 31. Plastics and their intermediates**
- 32. Pharmaceutical and biomedical applications**
 - Synthetic drugs
 - Pharmacokinetic studies
 - Drug monitoring
 - Toxicological applications
 - Plant extracts, herbal and traditional medicines
 - Clinico-chemical applications and profiling body fluids
- 33. Inorganic substances**
 - Cations
 - Anions
- 34. Radioactive and other isotopic compounds**
- 35. Other technical products and complex mixtures**
 - Surfactants
 - Antioxidants and preservatives
 - Various specific technical products
 - Complex mixtures and non-identified compounds
- 36. Thin-layer electrophoresis**
- 37. Environmental analysis**
 - General papers
 - Air pollution
 - Water pollution
 - Soil pollution
- 38. Chiral separations**

Cumulative CAMAG Bibliography Service (CCBS) Online Search

The screenshot shows the CCBS search interface. At the top, it says 'CUMULATIVE CAMAG BIBLIOGRAPHY SERVICE CCBS'. Below this is a search bar with the text 'Full text search :'. To the right of the search bar are two buttons: 'all editions' and 'search'. Below the search bar are three links: 'Classification', 'Keyword register', and 'CBS edition'. To the right of the search bar is a 'PDF Cart' icon with a shopping cart symbol and the text 'PDF Cart'. Below the search bar are five numbered callouts (1-5) explaining search options.

1 Full text search:
Enter a keyword, e.g. a substance name, a substance class, a technique, a reagent, or an author's name

2 Browse and search by CBS classification:
Select one of the 38 CBS classification categories and search by keyword

3 Alphabetical search:
select an initial character and browse associated keywords

4 Search by CBS edition:
Select a CBS edition and retrieve all abstracts published in this CBS issue. With this search you can get all abstracts of one CBS issue – similarly to the former printed yellow pages.

5 PDF Cart
Your cart is empty
Add items to your cart by clicking on the cart-icon to the left of the abstracts.
In the end, you can create a combined pdf document with all items put into the cart.

By using the cart icon you can create your individual selection of abstracts throughout CCBS search and export to PDF.

With the Cumulative CAMAG Bibliography Service (CCBS) Online Search, you can directly search for information within the CAMAG website. The CCBS covers more than 11'000 abstracts of TLC/HPTLC publications between 1982 and today. Providing reports was established by Dr. Dieter Jänchen in 1965 and the workflow and form evolved over the years. Manually written and transferred into the printed CBS form, it was set-up as a database in 1997. With CBS 93 in 2004, the electronic database for download was introduced and with the CBS 113 in 2014, the online search.

The database covers most relevant scientific journals in the field of Planar Chromatography including also various non-English publications in German, French, Spanish, Portuguese and Chinese. The CCBS features additional practical information for the analyst in the lab, for example details on the mobile phase or the detection. With CCBS the analyst is able to find relevant TLC/HPTLC publications which might be helpful for solving a particular analytical question.

Visit www.camag.com/ccbs and choose one of the following search options: full text search or search by CBS classification system or by alphabetical register or by CBS edition. For classical full text search, just enter a keyword in the search box, e.g., a substance name, a substance class, an analytical technique, a reagent, or an author's name, and find all related publications throughout the CCBS.


Alternatively, you can choose to browse by one of the 38 CBS classification categories and search by keyword. The alphabetical search allows selecting an initial character and browsing associated keywords. When browsing by CBS edition, you can retrieve all abstracts published in the corresponding CBS issue, formerly printed as yellow pages.

To create your individual selection of TLC/HPTLC abstracts add your preferred publications to the PDF cart and download these selected articles in one single PDF file.

CAMAG application tutorials



To illustrate the unique characteristics of high-performance thin-layer chromatography and to better support our customers in answering their analytical questions, we have recently started creating case studies that include tutorial videos for various HPTLC application fields.

These complementary videos visualize the practical steps of each HPTLC case study and present the required CAMAG instruments for performing the respective standardized method. All steps, including sample preparation, are shown. Go to www.camag.com/applicationfields, choose the preferred field of application and select a case study from the upper box on the right side. The following videos are already on our website:




Herbal medicines and botanical dietary supplements

- Quantification of ginkgolides A, B, C and bilobalide in *Ginkgo biloba* by HPTLC
- HPTLC-Fingerprint of *Ginkgo biloba* flavonoids



Food and feed stuff

- Quantitative determination of steviol glycosides
- Identification of fixed oils by HPTLC



Industrial applications

- Determination of monoacylglycerides in biodiesel



Forensics

- In-process control during chemical synthesis of ergoline psychedelics by HPTLC

Identification of herbal slimming drugs and screening for adulteration by HPTLC




From left: (back row) Kiran Dodtale, Shahrugh Bharucha, Dr. Sharad Medhe, Dr. Karuppiyah Jayachandran, (front row) Prashant Hande, Tukaram Thite

Anchrom has been CAMAG's distributor in India for more than 30 years. Mr. Thite runs Anchrom's well-equipped laboratory, established in Mumbai in 1989 with a team of experienced application specialists. His son, a medical doctor, came up with the idea to develop new HPTLC methods for identification and screening of herbal slimming drugs described here.

Introduction

Worldwide obesity has more than doubled since 1980. It substantially raises the risk for noncommunicable diseases such as cardiovascular diseases, diabetes, musculoskeletal disorders, and some cancers [1]. Medications available for weight loss are pharmacological agents that either reduce appetite or the resorption of food. In the United States one of the anti-obesity medications is orlistat, which has been approved by the FDA for long term use. It reduces intestinal fat absorption by inhibiting pancreatic lipase. Recently some herbal slimming drugs for weight reduction were commercially promoted, creating a concern in medical and scientific circles that there was a need to analyze samples of these drugs for synthetic drug adulteration. The draft USP chapter <2251> is specifically written for adulteration of dietary supplements (equivalent to herbal drugs) with drugs and drug analogs (for weight loss, sexual enhancement or sport performance enhancement).

HPTLC offers a fast way for identity testing of herbal slimming drugs and to check for adulteration. Herein, a combined methodology is proposed. The presence of orlistat can be determined prior derivatization by scanning densitometry and the identity can be confirmed by HPTLC fingerprint after derivatization.



Structure of orlistat

Chromatogram layer

HPTLC plates silica gel 60 F₂₅₄, 20 × 10 cm

Standard solutions

10 mg of orlistat dissolved with 10 mL acetone – methanol 2:3

Sample preparation

500 mg of each powdered herbal sample was mixed with 10 mL acetone – methanol 2:3 and sonicated for 15 min. After centrifugation the supernatant was used for application.

Sample application

Bandwise with Linomat 5, 15 tracks, band length 8 mm, track distance 11.4 mm, distance from left edge 20 mm, distance from lower edge 8 mm, application volume 10–20 µL for sample solutions and 5–10 µL for standard solutions

Chromatography

In the Automatic Developing Chamber 2 (ADC 2) with chamber saturation (with filter paper) for 20 min and after conditioning at 47% relative humidity using a saturated potassium thiocyanate solution, development with toluene – ethyl acetate 4:1 to the migration distance of 70 mm (from the lower edge), drying for 5 min.

Post-chromatographic derivatization

The plate was immersed into the phosphomolybdic acid reagent (10 g in 200 mL ethanol) using the Chromatogram Immersion Device, immersion time 3 s and immersion speed 3 cm/s and heating at 110 °C for 5 min using the TLC Plate Heater.

Documentation


With TLC Visualizer at 254 nm, 366 nm and white light illumination

Densitometry

TLC Scanner 4 and winCATS, absorption measurement at 195 nm, slit dimension 6.0 × 0.3 mm, scan speed 20 mm/s, spectra recording from 190 to 400 nm

Results and discussion

This method is well-suited for analysis of orlistat in herbal drugs, sold in tablet or powder form. The presence or absence of orlistat can be tested by scanning densitometry prior to derivatization.




Densitogram of different herbal samples (green and blue) and the orlistat standard (red, 10 µg/band)

To reach a low limit of detection (LOD), a high sample load can be applied on the HPTLC plate (e.g. 1 mg sample per 20 µL volume). As the LOD of the Orlistat standard is 70 ng/band, the LOD of the method is calculated to be 70 mg/kg. The normal dose prescribed for orlistat is 120 mg. Tablets usually weigh 1 g, while the recommended dose for powders is on average 4 g. Considering that about 12 % of orlistat were expected to be present in pills and about 3 % in powder samples, the method should be suitable for adulteration detection. In the case a zone is present at the same hR_F , a spec-

trum scan can be used to test for adulteration with orlistat. If illegally added Orlistat is found, HPTLC-MS can be used for further confirmation.

HPTLC is suitable for screening of a large number of samples. After derivatization, the HPTLC fingerprint can be used for further confirmation of orlistat and for comparing different herbal samples side by side.



HPTLC chromatogram under white light after derivatization with phosphomolybdic acid (samples U1–U5, orlistat standard S (10 µg/band) as white zone at hR_F 40)

Editor's note: The derivatization shown can be substituted by more appropriate reagents.

Anchrom and CAMAG will help interested parties to establish this method as a routine procedure for the screening of slimming drugs in herbal mixtures.

Further information is available on request from the authors.

[1] WHO, Fact sheet No. 311, Obesity and overweight, <http://www.who.int/mediacentre/factsheets/fs311/en/>

Contact: Dr. Sharad Medhe, Anchrom Testlab Pvt. Ltd., A-101, Shree Aniket Apt. Navghar Road, 400 081, Mulund (E) Mumbai, India, hptlc@anchrom.in

Determination of lactose in foodstuff



From left: Hans Griesinger, Michael Schulz, Markus Burholt, Vanessa Pilakowski, Michaela Oberle

TLC/HPTLC coupled with mass spectrometry has a wide range of application fields. Michael Schulz's group from Merck in Darmstadt develops new chromatographic layers for TLC and HPTLC. They also explore new fields of application for use with these and other techniques.

Introduction

Today food intolerances play an important role in our society, affecting the daily lives of many people. For example, a significant number of people are lactose intolerant, usually more in Asian countries than in Europe or North America. Foodstuff with a lactose content less than 100 mg/100 g can be labeled as "lactose-free", but currently no official limit has been established in Europe. Traditionally, analysis is done by HPLC, photometry or by an enzymatic test kit [1].

HPTLC offers an advantageous option for determination of lactose in high-matrix samples. A very easy and sensitive method was recently published [2]. As shown in this article, different food samples like milk and yoghurt can be applied on the HPTLC plate without a complex sample preparation. Lactose can be quantified by image evaluation with the help of the TLC

Visualizer and visionCATS. Afterwards the corresponding lactose zones can be eluted with the TLC-MS interface and detected by mass spectrometry.

Chromatogram layer

HPTLC plates silica gel 60 F₂₅₄ MS-grade (Merck), 20 × 10 cm

Standard solutions

Lactose in water (0.1 and 1 mg/mL)

Sample preparation

Milk was centrifuged for 5 min and the supernatant centrifuged again for 5 min. 1 mL was diluted in 50 mL water. 1 g yoghurt was dissolved in 10 mL water and centrifuged for 5 min. An instant sauce was cooked as directed. 1 g sauce was mixed with 10 mL water and centrifuged for 5 min.

Sample application

Bandwise with Automatic TLC Sampler (ATS 4), band length 6 mm, track distance 8.9 mm, distance from left edge 15 mm, distance from lower edge 10 mm, application volume between 0.3 and 2.0 µL/zone (samples) and between 0.1 and 4.0 µL/zone (standards)

Chromatography

In the Flat Bottom Chamber 20 × 10 cm with acetonitrile – water 3:1 (with 0.1 % trifluoroacetic acid) to the migration distance of 5 cm

Post-chromatographic derivatization

Spraying onto the HPTLC plate with aniline-di-phenylamine-phosphoric acid reagent (2 g diphenylamine and 2 mL aniline in 80 mL methanol, 10 mL of o-phosphoric acid 85 %, filled up to 100 mL) and heating on the TLC Plate Heater for 10 min at 120 °C

Editor's note: For quantification, a homogeneous application of the derivatization reagent is recommended. Automated immersion is better suited than manual spraying.

Documentation and digital evaluation


With the TLC Visualizer and *visionCATS* software under white light

Mass spectrometry

Elution of zones with TLC-MS Interface into an ESI-MS with acetonitrile – water 19:1 (with 0.1 % formic acid) and detected in positive ionization mode


Results and discussion

A representative densitogram of a yoghurt sample is shown. Lactose (2) shows the same hR_f values for the samples as for the standards and is well separated from matrix components (3).



Densitogram of yoghurt sample

The lactose content was determined by image evaluation with the help of the TLC Visualizer and *visionCATS*. Calibration was done in a range of 0.01 – 4 $\mu\text{g}/\text{zone}$.




Calibration curve: green standard levels, blue samples

| Determined lactose content | (g/100 g) |
|----------------------------|-----------|
| Organic fresh milk | 4.7 |
| UHT milk | 4.4 |
| Lactose free UHT milk | < 0.01 |
| Instant sauce | 1.6 |
| Yoghurt | 5.5 |

The left side of the plate was derivatized with aniline-diphenylamine-phosphoric acid reagent for localizing the zones on the right side of the plate (covered during spraying). The non-derivatized zones were eluted into the mass spectrometer.

Editor's note: It is an option one should consider to cut the developed plate with the smartCut and then derivatize one half of the plate.



Mass spectrum of the lactose zone in a yoghurt sample

The mass spectrum shows as mass signals of lactose (342 g/mol) the corresponding sodium adduct (m/z 365) as well as the solvent adduct (m/z 406). A semi-quantitative determination of lactose by image evaluation is fast and reliable. With this method lactose-free products can be confirmed as such.


Further information is available on request from the authors.

[1] L. Frank *et al.* GIT-Labor 15.01.2014

[2] G.E. Morlock, L.P. Morlock, C. Lemo, *J Chromatogr A* 1324 (2014) 215–223

Contact: Michael Schulz, Merck KGaA, Frankfurter Str. 250, 64293 Darmstadt, michael.schulz@merckgroup.com

Quantification of nicotine in liquids for electronic cigarettes



Tiên Do,
CAMAG

Introduction

The increasing availability of novel and emerging tobacco and related nicotine containing products, particularly electronic cigarettes (e-cigarette), will prompt authorities to create and implement new regulations. The legal status of e-cigarettes is currently pending in many countries. Some countries such as Singapore and Brazil have banned them. Because nicotine is an addictive and toxic substance, national legislations tend to be harmonized with respect to the levels of nicotine in the electronic liquid (e-liquid, the mixture used in vapor products including e-cigarettes). E-liquids are usually composed of nicotine in propylene glycol, glycerol, and water with added flavorings. In Europe consumers can buy e-liquids containing nicotine up to 20 mg/mL. Preliminary tests of some e-cigarettes conducted by the US FDA revealed inconsistencies of product labeling and nicotine content [1].

In anticipation of future legislation, HPTLC methods were developed for the quantification of nicotine by image evaluation and densitometry. The absence or presence of nicotine may be confirmed by hR_f values, UV and MS spectra compared to those obtained with the reference substance. To check compliance with the European limit of 20 mg/mL for nicotine in e-liquids, a limit test for controlling the maximum authorized content is proposed. Furthermore, the actual nicotine content can be determined based on a single-point or five-point calibration.

Chromatography layer

HPTLC plates silica gel 60 F₂₅₄ (Merck), 20 × 10 cm.

Sample preparation

Stock solutions of nicotine were prepared in propylene glycol/glycerol solution (1:1) at different levels from 0.135 mg/mL to 20 mg/mL. 0.3 g of each stock solution or sample is diluted with methanol to 10 mL.

Sample application

Bandwise with Automatic TLC Sampler (ATS 4), 15 tracks, band length 8 mm, distance from left edge 20 mm, distance from lower edge 8 mm, application volume 2 μ L

Chromatography

In the ADC 2 with chamber saturation (with filter paper) 20 min and after conditioning at 33 % relative humidity for 10 min using a saturated solution of magnesium chloride, development with toluene – acetone – diethylamine 10:10:1 to the migration distance of 70 mm (from the lower edge), drying for 5 min

Documentation

With the TLC Visualizer under UV 254 nm

Densitometry


TLC Scanner 4 and *visionCATS*, absorption measurement at 260 nm, slit dimension 5.0 × 0.2 mm, measurement speed 50 mm/s, spectra recording from 190 to 450 nm

Mass spectrometry


Elution of zones with TLC-MS Interface 2 (oval elution head) at a flow rate of 0.5 mL/min with methanol (with 0.1 % formic acid) into a single quadrupole mass spectrometer (ACQUITY QDa, Waters, USA) and detected in the positive ionization mode. Data processing and evaluation of mass spectra was performed with Empower.

Results and discussion

The developing solvent used for the separation of nicotine from other ingredients was modified [2]. The developing solvent was decreased in the elution strength. The hR_F value of 56 obtained for nicotine was sufficient for separation from other ingredients (propylene glycol, glycerol and flavors), which were located in the lower part of the chromatogram. Each sample was then visually evaluated for the presence of a zone matching the hR_F value of the nicotine. UV and MS spectra of target zones were recorded and compared with the one of the reference substance. Analysis of 9 samples including flavored and non-flavored e-liquids showed levels consistent with the respective label claims.



Left: Chromatogram under UV 254 nm (tracks 1–5: reference solutions, tracks 6–7: samples without flavor; tracks 8–9: samples with flavor); center: UV spectra of all reference and sample solutions; right: mass spectrum of the eluted nicotine zone in a sample at m/z 163 $[M+H]^+$



Calibration curve (image evaluation at UV 254 nm)

a) Limit test for samples at 20 mg/mL

A reference solution with 20 mg/mL nicotine (in propylene glycol/glycerol) and 9 e-liquid samples were applied. All samples show values below that of the reference solution.

| Sample | Claimed concentration (mg/mL) | Measured concentration (mg/mL) | |
|--------|-------------------------------|--------------------------------|---------|
| | | Scan* | Image** |
| 1 | 18.0 | 18.5 | 19.1 |
| 2 | 18.0 | 19.0 | 19.8 |
| 3 | 18.0 | 19.5 | 19.9 |

*extrapolation range of 10%

**extrapolation range of 15%

b) Single-point calibration

The consistency of the nicotine labeling of 3 samples, claimed to contain 18 mg/mL was tested. The measured concentrations were slightly higher, but below the European limit.

| Sample | Claimed concentration (mg/mL) | Measured concentration (mg/mL) | |
|--------|-------------------------------|--------------------------------|---------|
| | | Scan* | Image** |
| 1 | 18.0 | 18.5 | 19.1 |
| 2 | 18.0 | 19.0 | 19.8 |
| 3 | 18.0 | 19.5 | 19.9 |

*extrapolation range of 10%


**extrapolation range of 15%

c) Five-point calibration


For quantification of e-liquids containing different concentrations of nicotine a five-point calibration was used.

| Sample | Claimed concentration (mg/mL) | Measured concentration (mg/mL) | |
|--------|-------------------------------|--------------------------------|-------|
| | | Scan | Image |
| 1 | 6.0 | 6.5 | 5.7 |
| 2 | 12.0 | 12.3 | 12.6 |
| 3 | 18.0 | 17.7 | 17.8 |
| 4 | 18.0 | 17.2 | 17.9 |
| 5 | 18.0 | 17.9 | 18.1 |
| 6 | 18.0 | 16.6 | 16.7 |
| 7 | 18.0 | 16.4 | 16.2 |
| 8 | 18.0 | 16.5 | 16.4 |
| 9 | 18.0 | 16.6 | 17.6 |

The proposed method is applicable to e-liquid samples containing different concentrations of nicotine. No interference with excipients and flavor ingredients was found. Depending on the regulatory requirements to be met, analytical approaches can be adopted to either include image-based evaluation or densitometric determination of nicotine.



Single-point calibration curves of image-based or densitometric determination



Five-point calibration curves of image-based or densitometric determination

Further information is available on request from the author.

[1] Summary of Results: Laboratory Analysis of Electronic Cigarettes Conducted by the FDA: <http://www.fda.gov/NewsEvents/PublicHealthFocus/ucm173146.htm>

[2] Kenyon *et al.*, J AOAC Int 81 (1998) 44–50

Contact: Ti n Do, CAMAG, Sonnenmattstrasse 11, 4132 Muttenz, tien.do@camag.com

CAMAG TLC Scanner 4

Controlled by *visionCATS* HPTLC software, the CAMAG TLC Scanner 4 facilitates the densitometric evaluation of planar chromatograms.

The state-of-the-art software controls all functions of the TLC Scanner and enables optimal evaluation of measured data.




Technical specifications:

- Measurement of reflected light, either in absorbance or fluorescence mode
- Object formats up to 200 x 200 mm
- Spectral range 190–900 nm
- Automatic start of all lamps: deuterium, halogen-tungsten, and high pressure mercury lamp
- Data step resolution 25–200 μm
- Scanning speed 1–100 mm/s

visionCATS supports the following functions:

- Single-wavelength scan
- Multi-wavelength scan for quantitative evaluation of differently absorbing substances in one single measuring run
- Wavelength subtraction for background correction
- Spectra recording for checking identity
- Scanner selftest and IQ/OQ
- 21 CFR Part 11



Spectrum scan of sulfamerazine: overlay of sample and reference spectra

Further information:

www.camag.com/tlc-scanner

CAMAG[®]

WORLD LEADER IN PLANAR CHROMATOGRAPHY